

Version: 2.3.4

Running in align and count mode with the following parameters:

0 mismatch allowed

Minimal Phred Score per bp >= 30

Feature length: 20

Read alignment start position: 0

All data will be saved into C:/Users/afons/Downloads/bikard_fastq/2FAST2Q_output_2021_11_24_10_49_13\2FAST2Q_output_2021_12_20_13_17_45

Loading Features

11629 different features were provided.

Processing 118 files. Please hold.

Processing file 4 out of 118:	63% #####7		90870601/143945160	[00:05<00:03,	16673030.70characters/s]
Processing file 12 out of 118:	5% #1		14787364/298122739	[00:00<00:17,	16509475.69characters/s]
Processing file 6 out of 118:	50% #####3		87087627/175743839	[00:05<00:05,	16110963.00characters/s]
Processing file 11 out of 118:	9% ##2		26539623/284003755	[00:01<00:15,	16497925.05characters/s]
Processing file 3 out of 118:	79% #####7		89092403/112897070	[00:05<00:01,	16359410.16characters/s]
Processing file 5 out of 118:	58% #####5		90612488/155978997	[00:05<00:03,	16974129.56characters/s]
Processing file 9 out of 118:	37% #####2		86665482/234426059	[00:05<00:08,	16576674.86characters/s]
Processing file 7 out of 118:	47% #####6		88797884/189768453	[00:05<00:06,	16529244.99characters/s]
Processing file 10 out of 118:	35% #####3		90829008/262005192	[00:05<00:09,	17453406.44characters/s]
Processing file 8 out of 118:	44% #####		86629758/195837051	[00:05<00:06,	16187016.99characters/s]